



Phosphate, ortho, TNT, 0.06 - 5.00 mg/L PO₄

Product #:

2742545

Hazardous

AED Price:
Available

Contact Hach



Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests.

Specifications

Instrument:	DR 900
Measuring range:	0.06 - 5.00 mg/L PO ₄
Method:	Ascorbic Acid
Method number:	8048
Number of tests:	25 - 50
Parameter:	Reactive (Orthophosphate)
Platform:	Test 'N Tube

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)