



Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Product #:

212528

Hazardous

AED Price:

Contact Hach

No availability shown



Powder reagent in sealed packages for the determination of Ortho-Phosphate. For 25mL samples. pk/1000.

Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Measuring range:	0.02 - 2.50 mg/L PO ₄
Method:	Ascorbic Acid
Method number:	8048
Number of tests:	1000
Parameter:	Phosphate, ortho
Platform:	Powder Pillows
Quantity:	1000 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)