



PhosVer Phosphate Accuvac, 0.02-2.50 mg/L PO₄

Product #:

2508025

Hazardous

AED Price:
Available

Contact Hach



Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For determination of ortho Phosphate by the ascorbic acid method. Pack of 25.

Specifications

Instrument:	All except DR/820
Measuring range:	0.02 - 2.50 mg/L PO ₄
Method:	Ascorbic Acid
Method number:	8048
Number of tests:	25
Parameter:	Phosphate
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)