



Sodium Nitrate, 0.5M, 500 mL

Product #:	2507149
AED Price:	Contact Hach
Available	

Sodium Nitrate, 0.5 M, 500 mL, for 9245 and 9240 Sodium analysers

Specifications

Concentration:	5 mg/L
Instrument:	9245 & 9240
Parameter:	Sodium
Quantity:	500 mL
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)