



Hardness Indicator Solution, 5 - 100 mg/L, 500 mL

Product #:

2769249

Hazardous

AED Price:
Available

Contact Hach



For use with the SP510 Hardness analyser.

Specifications

Instrument:	SP510
Measuring range:	5 - 100 mg/L
Parameter:	Hardness
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)