



S5000 Phosphate Molybdovanadate Reagent (2.9 L)

Product #:

1420703

Hazardous

AED Price:
Available

Contact Hach



Specifications

Instrument: 5000 Phosphate

Measuring range: N/A

Parameter: Phosphate

Storage conditions: Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)