



Citric Acid / Surfactant Reagent for Silica analysis, 2.9 L

Product #:

2347003

Hazardous

AED Price:
Available

Contact Hach



Citric Acid / Surfactant Reagent (2.9 L) used in the Series 5000 Silica Analyser.

Specifications

Instrument:	5500sc Silica
Measuring range:	N/A
Parameter:	Silica
Quantity:	2.9 L
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)