



Anionic Surfactant Solution, 2.9 L

Product #:

2375503

Hazardous

AED Price:
Available

Contact Hach



For low range phosphate analyser, Series 5000.

Specifications

Instrument:	5000 Phosphate
Measuring range:	N/A
Parameter:	Phosphate
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)