



APA 6000 Hardness reagent 2, buffer solution, 1 L

Product #:

2695753

Hazardous

AED Price:
Available

Contact Hach



Specifications

Instrument:	APA6000
Measuring range:	N/A
Parameter:	Hardness
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)