



## APA6000 Low Range Hardness Standard 5 mg/L, 1 L

<b>Product #:</b>	<b>2696353</b>
AED Price:	Contact Hach
Available	

---

### Specifications

Concentration:	5 mg/L
Instrument:	APA6000
Parameter:	Hardness
Quantity:	1 L
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)