



APA6000 Alkalinity, Cleaning Solution, 1 L

Product #:	2697053
AED Price:	Contact Hach
Available	

The APA6000 Alkalinity Process Analyser uses five reagents each of which is supplied in a 1 L bottle designed to provide reagents for 30 days of continuous use.

APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample.

APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration.

APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration.

The fifth reagent, APA Alkalinity Cleaning Solution (269703), cleans the system during the Prime and Instrument Clean functions.

Specifications

Instrument:	APA6000
Measuring range:	N/A
Parameter:	Alkalinity
Quantity:	1 L
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)