



Total chlorine indicator solution for chlorine analyser CL17/CL17sc (473 mL)

Product #:

2263411

Hazardous

AED Price:
Available

Contact Hach



Specifications

Instrument:	CL17/CL17sc
Measuring range:	N/A
Parameter:	Chlorine, Total
Platform:	Solution
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)