



## Free chlorine buffer solution for chlorine analyser CL17/CL17sc (473 mL)

**Product #:** 2314111  
**AED Price:** Contact Hach  
**Available**

---

### Specifications

Instrument:	CL17/CL17sc
Measuring range:	N/A
Number of tests:	450
Parameter:	Chlorine, Free
Platform:	Solution
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)