



Free Chlorine Reagent Set for chlorine analyser CL17/CL17sc

Product #:

2556900

AED Price:
Available

Contact Hach

Hazardous



Formulated to meet the strictest quality standards

The Free Chlorine Reagent set is used with the CL17/CL17sc analyser as well as bulk analytical methods. Hach reagents provide a guarantee of accurate, consistent results.

Free Chlorine Reagent Set Consists of:

- One 24-gram bottle of DPD compound
- One 473-mL bottle of free chlorine indicator solution
- One 473-mL bottle of free chlorine buffer solution

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Specifications

Instrument:	CL17/CL17sc
Measuring range:	0.03 - 5.00 mg/L Cl_2 for CL17, 0.03 - 10.00 mg/L Cl_2 for CL17sc
Method:	DPD
Method number:	10059
Number of tests:	450
Parameter:	Chlorine, Free
Platform:	Solution
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
What's included?:	DPD compound, 24 g
	Free chlorine indicator solution, 473 mL
	Free chlorine buffer solution, 473 mL

What's included?

DPD compound, 24 gFree chlorine indicator solution, 473 mLFree chlorine buffer solution, 473 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)