



Stablcal Turbidity standard solution 1.0 NTU

Product #:

2659853

Hazardous

AED Price:
Available

Contact Hach



Stabilised Formazin turbidity standard for use in any turbidimeter. 1000 mL bottle.

Specifications

Concentration:	1.0 NTU
Measuring range:	N/A
Parameter:	Turbidity
Quantity:	1 L
Storage conditions:	Store between 5 - 25 °C, away from direct sunlight

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)