



Stablcal Turbidity standard solution 20 NTU (1000 mL)

Product #:

2660153

Hazardous

AED Price:
Available

Contact Hach



Stablcal Turbidity Standard for all Turbidimeters. USEPA accepted for calibration purposes. Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. NTU Value: 20.0. Bottle: 1000 mL (1 L).

Specifications

Concentration:	20 NTU
Measuring range:	N/A
Parameter:	Turbidity
Quantity:	1 L
Storage conditions:	Store between 5 - 25 °C, away from direct sunlight

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)