



Turbidity standard 4,000 FNU, 100 mL

Product #:

LCW813

Hazardous

AED Price:

Contact Hach



Specifications

Concentration:	4000 FNU
Measuring range:	N/A
Parameter:	Turbidity
Storage conditions:	Store in a cool, dark place, away from direct sunlight

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)