



Sension+ ORP/Redox Gel electrode, ceramic diaphragm, for portable meter

Product #: LZW5055.97.0002
AED Price: Contact Hach
Available

ORP/Redox electrode with platinum wire. Polycarbonate body. Range ± 2000 mV, 0-80 °C. Low maintenance. For portable meters.

Specifications

Accuracy:	1 mV with pH1 Meter
Cable length:	1 m
Dimensions (D x L):	12 mm x 85 mm (w/o handle)
Electrode junction :	Ceramic pin
Electrode type:	Portable
Filling solution:	Non-refillable gel
Length:	85 mm (w/o handle)
Material:	Sensor Body: Polycarbonate
Measuring range:	± 2000 mV
Parameter:	ORP
Sensor type :	Platinum annular ring
Temperature range:	0 - 80 °C
Type:	Portable
Warranty:	6 months