



CDC745-9 General-purpose conductivity cell, 2-pole

| | |
|--------------------|--------------|
| Product #: | E61M013 |
| AED Price: | Contact Hach |
| Call for ship date | |

2-pole conductivity electrode with screw cap and epoxy body. Cell constant 1cm^{-1} . Platinised, without temperature sensor. Max. temperature 100°C , min. immersion depth 14mm.

Specifications

| | |
|--------------------|---|
| Conductivity Cell: | CDC745-9 |
| Electrode: | 2-pole |
| Special Features: | Conductivity electrode with screw cap and epoxy body. |
| Warranty: | 12 months |