



Boron Reagent Powder Pillows, 0.2-14.0mg/L B

Product #:	1417099
AED Price:	Contact Hach
Available	

Powder reagent in sealed packages for the determination of Boron with Carmine Method. pk/100.

Specifications

Footnote:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	All except DR/800, PC II
Measuring range:	0.2 - 14.0 mg/L B
Method:	Carmine
Method number:	8015
Number of tests:	100
Parameter:	Boron
Quantity:	100 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)