



Molybdate Reagent Powder Pillows

Product #:	104199
AED Price:	Contact Hach
Available	

Powder reagent in sealed package for use in the Silicomolybdate Method for the determination of high range Silica (0-100.0mg/L SiO_2). For 25mL sample. Pack of 100.

Specifications

Method:	Silicomolybdate
Method number:	8185
Number of tests:	100
Parameter:	Silica
Platform:	Powder Pillows
Quantity:	100 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)