



Chromium Reagent, Accuvac, 0.010-0.700 mg/L Cr (VI)

Product #: 2505025

Hazardous

OBSOLETE ITEM
This item is no longer available.



Suggested replacements:

- 1271099

Ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Chromium, hexavalent. Pack of 25 ampules.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/850, DR/820
Measuring range:	0.010 - 0.700 mg/L Cr ⁶⁺
Method:	1,5-Diphenylcarbohydrazide
Method number:	8023
Number of tests:	25
Parameter:	Chromium
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)