



TitraVer EDTA hardness titrant, 0.020 N, 1L

Product #:	20553
AED Price:	Contact Hach
Available	

Titrant for the determination of hardness, 0-5000 mg/L CaCO₃. 1 Liter bottle.

Specifications

EPA compliant:	Yes
Measuring range:	0 - 5000 mg/L CaCO ₃
Method:	EDTA Titration
Number of tests:	40 - 100
Parameter:	Calcium and Total Hardness
Platform:	Solution
Quantity:	1000 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)