



BOD nutrient buffer pillows, 0.5 mL size for 300 mL dilution water

Product #: 1416066
AED Price: Contact Hach
Call for ship date

BOD Reagents

For preparation of BOD5 dilution water. Each pillow combines the APHA-specified buffer and nutrients to prepare 300 mL of dilution water. Pack of 50 pillows.

Specifications

Application:	Dilution water preparation
Description:	Nutrients & Buffer
Formula:	Nutrients & Buffer
Package Type:	Powder Pillows
Parameter:	BOD
Platform:	Powder Pillows
Quantity:	300 mL
Units:	50/pk

Required Accessories

- BOD Seed Inoculum, Polyseed, pk/50 capsules (Item 2918700)
- BOD Standard Solution, 300 mg/L, pk/16 - 10-mL Voluette Ampules (Item 1486510)
- Nitrification inhibitor for BOD (35 g) (Item 253335)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)