



Filling Solution, Reference, 3 M KCl with AgCl, 100 mL

Product #:	S21M004
AED Price:	Contact Hach
Ships within 2 weeks	

Filling solution for electrode reference element, Potassium Chloride (KCl) with Saturated AgCl, 3 M, 100 mL bottle

Specifications

Concentration:	3 M KCl, saturated with AgCl
Description:	MTC301, PHC301
Formula:	KCl 3M & AgCl Sat.
Parameter:	pH, ORP
Platform:	Solution
Quantity:	100 mL
Solution Type:	Refill