



APA6000 Low Range Hardness Standard, 50 µg/L, 1 L

Product #: 2696253
AED Price: Contact Hach
Available

Specifications

Concentration:	50 µg/L
Instrument:	APA6000
Parameter:	Hardness
Quantity:	1 L
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)