



## M291Ag Metal electrode for silver halides titration

<b>Product #:</b>	<b>E34M002</b>
AED Price:	Contact Hach
Ships within 2 weeks	

With screw cap. Sensing element silver rod, salt-bridge solution saturated  $K_2SO_4$ . Temperature range 0-80°C.

---

### Specifications

Electrode type:	M291Ag
Warranty:	6 months
What's included?:	Includes: Radiometer Analytical M291Ag-9 metal electrode with Certificate of Conformity.

---

### What's included?

Includes: Radiometer Analytical M291Ag-9 metal electrode with Certificate of Conformity.