



## Nitrification inhibitor for BOD, formula 2533(TM), TCMP, 500 g

**Product #:** 253334  
**AED Price:** Contact Hach  
**Available**

### BOD Reagents

TCMP - 2-chloro-6 (trichloromethyl) pyridine (N-Serve), coated on an inert substrate. Use 0.16 gr in a 300 mL sample to inhibit nitrogenous oxygen demand in the BOD test 500 g.

Alternatively, use with bottle PN 1704200 and dispenser cap, Part No. 45901.

---

### Specifications

Application:	Nitrification Inhibitor (CBOD)
Description:	Nitrification Inhibitor
Formula:	Trichloromethyl Pyridine (TCMP)
Method number:	8043
Number of tests:	2500
Package Type:	Bottle
Platform:	Jar
Quantity:	500

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)