



ez GGA - No Pipette GGA, Single Dose BOD Standard, pk/20 2-mL ampules

Product #: 2514420
AED Price: Contact Hach
Available

BOD Reagents

ez GGA, conveniently packaged in single-use ampules. No more diluting, no more pipetting, and no more guessing. Simply crack the top of the ampule open and drop the entire ampule into your BOD bottle. The pre-measured GGA contains the standard 198 mg/L of BOD from Glucose plus Glutamic Acid (GGA). The presence of the glass ampule in your BOD bottle serves as a quick visual indicator that you added the GGA to your sample. Product is equivalent to Standard Methods.

Specifications

Application:	Calibration Standard
Concentration:	150 mg/L
Description:	ez GGA pre-measured Glucose plus Glutamic Acid (Dilution Method)
Formula:	Glucose & Glutamic Acid (GGA)
Package Type:	Ampuls
Platform:	Solution
Quantity:	mL/20 ez GGA Ampules