



Filling Solution, Internal, Polarographic DO, 30 mL

Product #: 2759123
AED Price: Contact Hach

Filling solution for Polarographic Dissolved Oxygen probe, Potassium Chloride/Glycerine, 30 mL bottle

Specifications

Description:	5130, 5197000
Formula:	KCl & Glycerine
Manual languages:	Okay. 16.10.2012 Bauer
Quantity:	30 mL
Solution Type:	Refill

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)