



## Pipet tips, TenSette pipet 1970001, 0.1-1.0 mL, pk/1000

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2185628</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Perform accurate, economical pipetting with Hach's Tensette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

---

### Specifications

|           |          |
|-----------|----------|
| Quantity: | 1000 /pk |
|-----------|----------|