



## M241Pt Metal Electrode, Platinum Plate, d=7.5 mm, Banana Plug (Radiometer Analytical)

<b>Product #:</b>	<b>E31M001</b>
AED Price:	Contact Hach
Available	

The Radiometer M241Pt Metal Electrode has a platinum plate sensor and is recommended for general purpose measurements. An external reference electrode is required. The M241Pt has a 1 meter fixed cable with banana plug.

### The right electrode for your application.

To meet your analysis needs, the Radiometer Analytical electrode range features single, double, and combined metal electrodes with platinum, silver, gold, antimony, mercury and glassy carbon sensing elements. This variety of designs can be used for redox measurements, redox titrations, potentiometric techniques, and imposed current potentiometric titration.

---

## Specifications

Electrode type:	M241Pt
Warranty:	6 months
What's included?:	Includes: Radiometer Analytical M241Pt metal electrode with Certificate of Conformity.

---

## What's included?

Includes: Radiometer Analytical M241Pt metal electrode with Certificate of Conformity.