



pHpulp PT100 pH Glass Electrode with 5 m Fixed Cable

Product #: LZX475
AED Price: Contact Hach

With integrated temperature sensor (PT100).

Specifications

Cable length:	5 m fixed cable
Diaphragm:	Ring diaphragm, porous PTFE, double
Electrolyte:	Exterior EPH gel, interior standard gel KCl/AgCl
Length:	300 mm
Material:	Stainless steel
Measuring range:	0 - 14 pH
Operating temperature range:	-5 - 40 °C
Pressure range:	10 bar at 25 °C
Reference system:	Ag/AgCl
Special notes:	Suitable for use with the 22 mm inline armature (LZX 467)
Temperature sensor:	PT100
Warranty:	6 months