



## LTLCONT pH Electrode with Top68 Connector

**Product #:** LZX536  
**AED Price:** Contact Hach

Low temperature glass and electrolyte, single reference system, electrode for ion-weak and/or cold water (10-100 $\mu$ S/cm), without temperature sensor

---

### Specifications

Cable connection:	TOP68 connector
Diaphragm:	Ring diaphragm, porous PTFE
Electrolyte:	KCl, saturated KCl/AgCl (crystal)
Immersion depth:	105 mm
Length:	140 mm
Measuring range:	2 - 11 pH
Operating temperature range:	-20 - 50 °C
Pressure range:	3.5 bar
Reference system:	Ag/AgCl
Sensor thread:	PG 13.5
Temperature sensor:	None
Warranty:	6 months