



SOTA PT100T pH electrode, 12 mm, TOP68 connector

Product #: LZX583
AED Price: Contact Hach

Standard glass and reference system with head connector TOP68

Specifications

Application:	General use
Cable connection:	TOP68 connector
Connector:	TOP68
Diaphragm:	Ring diaphragm, porous PTFE, double
Electrolyte:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Immersion depth:	105 mm
Length:	140 mm
Material:	Sensor: Glass
Measuring range:	0 - 80 °C
Pressure range:	0 - 3.5 bar
Reference electrode:	None
Sensor thread:	Pg 13.5
Temperature sensor:	Pt100
Warranty:	6 months