



pH sensor without integrated temperature sensor

Product #: LZY023
AED Price: Contact Hach

pH combination electrode with PG 13,5 screw plug connect mounting, polymer electrolyte and non refiliable KPG diaphragm; Pressure Equipment Directive for Gas 1/Liquids 1 Art. 3.3 Sep

Specifications

Cable connection:	S7 connector
Diaphragm:	Annular gap
Electrolyte:	Polymer with visible KCl reserve
Length:	120 mm
Measuring range:	2 - 13 pH
Operating temperature range:	0 - 100 °C
Pressure range:	10 bar
Reference electrode:	Ag/AgCl
Reference system:	Ag/AgCl
Response time:	< 30 s
Temperature sensor:	None
Warranty:	6 months
Weight:	0.15 kg