



Analogue combination ORP sensor, 3/4" NPT, 10 m cable

Product #: 08351=C=0000

AED Price: Contact Hach

Material PPS, Platinum ORP electrode

Specifications

Application:	ORP for high temps
Cable length:	10 m fixed cable
Connector:	Fixed cable
Drift:	<3 mV/week
Electrolyte:	Polymer electrolyte, KNO ₃ and KCl
Impedance:	Reference impedance @ 25 °C: <50 kΩ
	Glass impedance @ 25 °C: Platinum
Length:	150 mm
Material:	Reference Electrode: Ag/AgCl
Measuring range:	-1500 mV - 1500 mV
Operating temperature range:	Temperature Maximum: 110 °C (230 °F)
Pressure range:	Maximum Pressure: 10 bar (145 psi) at 80 °C
Process connection:	0.75 " NPT thread at both ends
Reference electrode:	None
Slope Range:	N/A
Temperature sensor:	PT100
Warranty:	6 months
Weight:	0.62 kg