



Polymetron 8429B pH reference electrode

Product #: 368429,00000
AED Price: Contact Hach

Single reference with KCl gel liquid junction is a ceramic frit, 0 - 80°C, R < 10 KΩ, to be mounted on an electrolytic bridge 8429.1 (PP) or 8429.2 (PVDF)

Specifications

Cable connection:	S7 connector
Length:	30 mm
Load:	< 10 KΩ
Operating temperature range:	0 - 80 °C
Special notes:	single reference with KCl gel; liquid junction is a ceramic frit