



Polymetron 8444B Standard reference electrode (8444B)

Product #: 368444,00000

AED Price: Contact Hach

Specifications

Cable connection:	S7 connector
Diameter:	12 mm
Diaphragm:	Ceramic frit
Immersion depth:	117 mm
Impedance:	< 20 k Ω
Operating temperature range:	-5 - 120 °C
Process connection:	tubing for electrolyte supply (tubing not included)
Sensor thread:	Pg 13.5
Special notes:	Refillable electrode, permits pressurized electrolyte
Warranty:	6 months