



## Replacement pH electrode for digital 1200-S sc pH sensor

**Product #:** LZX889  
**AED Price:** Contact Hach

---

### Specifications

Accuracy:	$\pm 0.02$ pH
Cable connection:	S7 plub connector
Diameter:	12 mm
Material:	Glass (general purpose)
Repeatability:	$\pm 0.5$ %
Response time:	< 15 s
Sensitivity:	$\pm 0.5$ %
Special notes:	Only applicable in combination with 1200-S sc pH sensor
Warranty:	6 months