



Pole mounting hardware pH, 10 cm bracket, stainless steel pole 2 m

Product #: LZY714.99.23920

AED Price: Contact Hach

Stainless steel pole mounting hardware for 1" NPT pH sensor including 10 cm stainless steel bracket and 2 m stainless steel pole. Suitable for all 1" NPT pH sensors like 1200-S sc (LXV426.99.10001) and pH-D-S (LXV427.99.10001).

Specifications

Diameter:	40 mm
Instrument:	Suitable for all 1" NPT pH sensors
Length:	Pole: 2000 mm
Material:	Wetted parts: pole SS 1.4301; sensor adapter PVC; screw set SS.
Weight:	Lifting part only: 4.9 kg