



3M KCl electrolyte for reference electrode, 500mL

Product #:	Z363140,00500
AED Price:	Contact Hach
Call for ship date	

Consumables for a sodium analyser in its spare kit.

Specifications

Concentration:	3.0 mol/L
Quantity:	500 mL
Special notes:	Reference electrolyte for refillable reference electrodes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)