



Ammonium/Potassium Combination Standard Solution: 5.000 mg/L NH₄-N, 20.000 mg/L K, 1000 mL

Product #: LCW892
AED Price: Contact Hach

Ammonium/Potassium standard for validation purpose of ammonium process probes.

Specifications

Concentration: 5000 mg/L NH₄-N, 20.000 mg/L K
Instrument: Amtax sc
Parameter: Ammonium

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)