



Amtax sc CAL1: Standard solution 50 mg/L NH₄-N (1.9 L), MR3: 10-1000 mg/L NH₄-N

Product #: BCF1012
AED Price: Contact Hach

Specifications

Concentration:	50 mg/L
Instrument:	Amtax sc
Measuring range:	10-1000 mg/L NH ₄ -N
Parameter:	Ammonium

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)