



Amtax inter2 Standard solution 0.5 mg/L $\text{NH}_4\text{-N}$, 1000 mL

Product #:	LCW862
AED Price:	Contact Hach

For MR1: 0.02 - 2.00 mg/L $\text{NH}_4\text{-N}$

Specifications

Concentration:	0.5 mg/L
Instrument:	Amtax inter 2
Parameter:	Ammonium

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)