



Reagent C for LCW893, Phosphax Sigma

Product #:

BCF691

Hazardous

AED Price:

Contact Hach



Specifications

Instrument:	Phosphax sigma
Measuring range:	0.01 - 5 mg/L TP
Parameter:	Total Phosphorous

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)