



Sigmatrax 2 Sampling Unit, with 20 m connecting hose, w/o control unit

Product #: LXV232
AED Price: Contact Hach

Using the immersed sampling probe the volume of sample necessary for the analysis is pumped without coming into contact with pumps.

The sample is delivered to a small glass vessel in the control unit and is homogenized using an ultrasonic generator.

Specifications

Delivery height:	max. 7 m
Diameter:	184 mm probe
Hose length:	20 m
Length:	404 mm probe
Operating temperature range:	5 - 30 °C Sample
Power consumption:	350 VA
Warranty:	24 months
Weight:	15 kg probe with 20m hose
What's included?:	Sampling Unit only, requires Sigmatrax 2 control unit

What's included?

Sampling Unit only, requires Sigmatrax 2 control unit

Required Accessories

- Sigmatrax 2 Control unit sampling & ultrasonic homogenisation (Item LXV215)