



Iron trace cuvette test 0.01-1.0 mg/L Fe, 20 tests

Product #:

LCK521

Hazardous

AED Price:

Contact Hach



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Convenient handling with 50 mm rectangular cuvettes + Dosicaps

Highest accuracy for analysis of ultra low range concentrations

Specifications

According to standard:	ISO 6332-1988, DIN 38406 E1-1
Description:	Iron trace 50 mm cuvette test
Digestion required:	Yes, for Total Iron
Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.01 - 1.0 mg/L Fe
Method:	1.10-Phenanthroline
Number of tests:	20
Parameter:	Iron, trace
Platform:	LCK
Special remarks:	For total Iron digestion use Crack Set LCW902
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)