



Manganese Reagent Set 0.005-0.5 mg/L Mn

Product #:

LCW532

AED Price:

Contact Hach

Hazardous



For determination of Manganese by the PAN method.

Specifications

Measuring range:	0.005 - 0.5 mg/L Mn
Method:	1-(2-pyridylazo)-2-naphthol (PAN)
Number of tests:	50
Parameter:	Manganese
Special remarks:	50mm rectangle-cuvette L郑341 needed
Storage conditions:	15 °C - 25 °C

Required Accessories

- Sample cell, 50 mm rectangular PMMA, 10 pcs. (Item L郑341)