



Crack Set Reagent set for metal digestions

Product #:

LCW902

AED Price:

Contact Hach

Hazardous



For digestions of complex, colloidal or undissolved metals in water and waste water

Specifications

Description:	Sample preparation for the determination of total heavy metals with LCK306 Lead, LCK308 Cadmium, LCK321 Iron, LCK329 Copper, LCK337 Nickel, LCK360 Zinc, LCK521 Iron trace, LCK529 Copper trace, LCK537 Nickel trace
Method:	Acid digestion with oxidising agent
Number of tests:	50
Parameter:	For Lead, Cadmium, Iron, Copper, Nickel, Zinc
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)